

SALBUTAMOL + GUAIFENESIN

PLEURAVENTIN 1 mg / 50 mg per 5 mL SYRUP ANTI-ASTHMA / EXPECTORANT

PHARMACOLOGY :

Salbutamol is a selective beta 2- adrenoceptor agonist. At therapeutic doses it acts on the beta 2-adrenoceptor of bronchial muscles with or no action on the beta 1-adrenoceptors of the heart. Guaifenesin can make the viscous mucus of the respiratory pathway more fluid and therefore facilitates expectoration and reduces cough. In the presence of excessive mucus adequate pulmonary ventilation cannot be achieved simply by relieving bronchospasm, it is necessary to remove mucus obstruction. The combination of salbutamol and guaifenesin is designed to relieve respiratory obstruction and improve pulmonary ventilation.

FORMULATION:

Each 5 mL (one teaspoonful) contains:
Salbutamol 1 mg
(equivalent to Salbutamol sulfate 1.2 mg)
Guaifenesin 50 mg

INDICATION:

Productive cough as in bronchitis, bronchial asthma, pneumonia and other acute respiratory tract infections.

CONTRAINDICATION:

PLEURAVENTIN Syrup is contraindicated in patients with a history of hypersensitivity to any of its component.

ADVERSE EFFECTS:

Headache, urticaria, hypotension, feeling of tension, hand tremor, Peripheral vasodilation, transient and muscle cramps.

DOSAGE:

Two or three times a day.
Adults : 2 to 4 teaspoonfuls
Children:
6 to 12 years old : 2 teaspoonfuls
Under 6 years old : 1 to 2 teaspoonfuls. Or as prescribed by the physician.

CAUTION:

Foods, Drugs, Devices and Cosmetics Act prohibits dispensing without prescription.

STORAGE CONDITION:

Store at temperatures not exceeding 30°C

AVAILABILITY:

Amber bottle /Alu cap plain of 60 mL and 120 mL

SHELF LIFE : 36 Months

Manufactured by
SAN MARINO Laboratories CORP
1 Crisanto delos Reyes Street, Brgy. Javalera
General Trias, Cavite
For:
JOHNTANN INTERNATIONALE PHARMA CORP.
25 Kabignayan Street, Quezon City